

whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support whenever you need support

product catalogue



whenever
you need
support ...



İZMİR ORTOPEDİ

Orthopaedic Soft Goods, Inc.



Izmir Ortopedi was founded in 1970 by Altan Penol and two workers as a small workshop. Increasing the number of its new products every year, with its registered brands like Altex, Erüs and Izor, the company is producing over 200 items of products in the year 2003 with its 110 workers. Izmir Ortopedi, as being the biggest orthopedic soft goods producer in Turkey, is setting the standards of the Market.

With this responsibility in mind, by registering its standards, Izmir Ortopedi with its workers distributors and suppliers is increasing its efforts day by day for being a good player in the world league, and besides with ISO 9001 standards the company has targeted to maximize the customer satisfaction by applying a management system that disciplines all the related processes beginning with the research of the customer demands and the market and leading to design, production and after sales service.



Marketing

With its active sales department Izmir Ortopedi has a face to face contact with customers, delivers instantly to company addresses and serves all over the country.



Communication

Our primary importance is servicing. On the other hand of the phone line team of experts is always ready to inform you about our product groups and innovations. We will also be happy to meet you face to face at your address on your demand.



On-line

In endless possibilities of communication we are always here whenever you wish to reach us. With our continously updated web site be informed of our innovations or order on-line if you wish.

www.izor.com.tr

- ➔ 3400 m2 of indoor production facility.
- ➔ A staff of 110 specially qualified on their own branches.
- ➔ A widespread service and distributor net among Turkey
- ➔ Over 200 items product range
- ➔ New investments everyday for the development of the market.

ISO CERTIFIED MANUFACTURER



SINCE 1970



capacity



technology



quality



OEM DESIGNS ARE WELCOME

Erüs® *Medical Corset*

Foot Care
Hand-Wrist Orthoses
Corsets
Cervical Collars-Braces
Various Products
Hip Orthoses
Orthopedic Shoes
Orthotics-Insoles
Silicone Products

6-7
8-12
13-18
19-22
23
24
25
26
27-29

Neoprene Supports **ALTEX®**

Foot-Leg
Waist-Back
Hand Arm-Shoulder
Hot-Cold Compress

32-34
35-37
38-40
41-45
46

ALTEX® *Medical Stockings*

Medical Stockings

48-49

iz/DR® Health Care Products

Walking Aids
Wheelchairs

52-53
54-57


Erüs Size Chart
Altex Size Chart

58-59
60-61

medical corsets

every moment of life

ERÜS[®] *Medical Corset*



❖ METATARSAL PAD

301-001



- Placed proximal to metatarsal heads to unload pressure across metatarsal arch and relieve metatarsalgia. constructed with the pad placed behind the ball-of-the-foot to relieve pressure and redistribute weight from the painful area to more tolerant areas.

■ Indication
Metatarsalgia



❖ METATARSAL PAD

301-186



- This model is more soft than 301-001 and fits every foot. It's fitted to second toe. Designed for Adults.

■ Indication
Metatarsalgia



❖ ANKLE SUPPORT

301-006



- Helps protect ankles from further trauma. Wraps under arch and around ankle for non-slip placement. Stretched to fit most adults. 8 size available.

■ Indication
Ankle sprains



❖ ANKLE SUPPORT

301-007



- Support for weak or injured ankles. Comfortable to wear with a contoured Fit. This support is designed to supplement the supporting action of the ankle's muscles and ligaments while they regain strength after injury.

■ Indication
Ankle sprains



301-077

ANKLE SPLINT

- This lightweight ankle splint may be used for either therapeutic or preventative support for the ankle. It is effective in stabilizing the ankle joint following severe sprains, tendinitis or after cast removal. The rigid pre-bent metal stays are adjustable and provide excellent side-to-side support. The stays also help the splint control rotational strains on the ankle, and provide resistance to inversion and eversion.

Indication

Tendon injuries, Ankle sprains



304-021

HALLUX VALGUS DAY BANDAGE

- Hallux Valgus Day Bandage were designed to relieve the pain caused by bunions. The unique sling design relieves discomfort by gently realigning the big toe and protecting the bunion area with a soft cushion. Simply slip the elastic over the instep of the foot and place the nylon toe cap over your big toe. Then pull the elastic strap around the back of your heel and fasten. Can be used daytime with shoes. There is a metatarsal pad bottomsides.

Indication

Hallux Valgus



304-035

HALLUX VALGUS NIGHT SPLINT

- The Hallux Valgus Night Splint helps keep your big toe aligned to ease pain caused by bunions. The soft bumper pad cushions the big toe while the velcro strap acts as a lever to gradually align the big toe relieving discomfort. Simply slip the toe piece between toes and slide sleeve over bunion. Tightening the velcro strap then open the space between your first and second toe. One size.

Indication

Hallux Valgus



304-129

HALLUX VALGUS TOE SPREADER

- Relieves between toe friction and irritation. Maintains toe alignment post digital or forefoot surgery. Soft silicone. Small and large size available.

Indication

Hallux Valgus



ALUMINIUM FINGER SPLINTS



- 1 Model 500 Finger Splint
- 2 Model 502 T Splint
- 3 Model 502-001 Baseball Splint
- 4 Model 502-002 Curved Finger Splint
- 5 Model 502-003 Finger Cot
- 6 Model 502-004 Finger Protector
- 7 Model 502-005 Frog Splint

■ Indication

Finger Injuries, Finger Fractures, Muscular or tendon contractions.
Dislocations / Sprains of the PIP and DIP joints.
For treating scars



MALLET FINGER SPLINT

508



- Moulded polyethylene splints to support the distal joint in extension. Secure with strapping.

■ Indication

Finger Injuries, Finger Fractures, Mallet Finger



WRIST SUPPORT

504



- Help maintain normal hand activity after sprains and with tendinitis. Used to help prevent strain and injury during work or athletics. Easy off - Easy on velcro closure.

■ Indication

Sprains, Tendinitis



THUMB SPLINT

504-001



- Helps support and limit the motion of the MP joint of the thumb. Ideal for use on soft tissue injuries, ligament strains, Gamekeeper's thumb, osteoarthritis, and degenerative joint disease. Allows for full finger dexterity. Thermal lining provides warmth to the thumb area and joint..

■ Indication

Gamekeeper's Thumb, CMC and MC injuries



505

WRIST SPLINT

- It's used when the wrist is strained and twisted or when the wrist needs to have a rest. Designed to stabilize the wrist in the neutral position. This splint stabilizes the wrist while still allowing full movement of the fingers. Features a removable palmar stay for extra support and immobilization.

Indication

Rheumatic and traumatic arthritis, tendinitis and carpal tunnel, post op



505-001

WRIST SPLINT WITH ABDUCTED THUMB

- Recommended when immobilization of the wrist and/or thumb is needed. Provides intimate support, immobilization and thumb abduction.

Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, Gamekeeper's Thumb, CMC and MC injuries



506

FOREARM WRIST SPLINT

- The splint's longer forearm length comfortably distributes pressure and weight over two-thirds of the forearm. Wrist motion is restrained without interfering with critical joint motion.

Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, CMC and MC injuries, when wrist and forearm needs to have rest



516

TENNIS ELBOW

- This armband applies focused compression to the tendon without restricting circulation. Eliminates stress and helps reduce inflammation and pain. Ideal for wear during athletic activity, work, or at rest. One size.

Indication

Epicondylitis



EPICONYLITIS BANDAGE

517



- Helps relieve pain from muscular overuse caused by work, sports or injury. Provides targeted compression to the long tendons of the forearm allowing them to remain parallel and pain free. Dual action compression works simultaneously to treat lateral and medial epicondylitis (forehand and backhand tennis elbow). Ideal for epicondylitis, tennis elbow and muscle strains.

■ Indication

Carpal Tunnel Syndrome, arthritis, and tendonitis



THUMB SPICA

520



- Restricts thumb motion and supports CMC and MP joint arthritis. Allows gripping action for ski pole, hammer handle, etc. Fabricated from polyethylene in a unique and comfortable design. A single soft hook and loop strap closure allows for ease of application. Left or Right S-M-L

■ Indication

Gamekeeper's Thumb, CMC and MC injuries



DE QUERVAIN SPLINT

521



- Anatomically shaped of low density polyethylene with closed-cell foam lining. Immobilizes and protects injured thumb joints.

■ Indication

For sprains of thumb, soft tissue injuries, gamekeeper injuries, collateral ligament support, and avulsion fractures of the thumb



WRIST IMMOBILIZATION SPLINT

522



- Anatomically shaped of low density polyethylene with closed-cell foam lining. Ribbed tongue rigid yet comfortable for most traumatic, post-surgical, or longer term rehabilitation applications. Adjustable velcro straps allows for swelling and facilitate wound management.

■ Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, CMC and MC injuries, when wrist and forearm needs to have rest



523

FUNCTIONAL WRIST SPLINT

- Keeps the functional position of the hand. Fabricated from polyethylene in a unique and comfortable design

Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, CMC and MC injuries, when wrist and forearm needs to have rest



524

POSITIONAL HAND SPLINT

- Keeps the wrist in neutral position, thumb in abduction position. Device for progressive and total extension of the hand.

Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, CMC and MC injuries, when wrist and forearm needs to have rest



510 - 512

ARM SLING

- Model 510: made of soft, breathable cotton fabric that is light weight and soft to the touch. Slide buckle on strap for easy adjustment .
- Model 512: waist belt for immobilization included.

Indication

Sprained, broken or surgically operated arm while it recuperates



513

FOAM ARM SLING

- Superior arm sling that evenly distributes weight, while holding arm securely in place. Made of durable foam padded cotton with hook and loop closure for easy adjustment and maximum comfort.

Indication

Sprained, broken or surgically operated arm while it recuperates



ARMSLING 300 ABDUCTION

514



- Can be used after rotator cuff repair. Helps keep the shoulder in slight 30° abduction. Used as a transition between immobilization with a shoulder abduction splint and full function without a splint.

■ Indication

For bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



SLING AND SWATHE IMMOBILIZER

518



- This sling gives lightweight protection in a cradletype style that fits either right or left arms. Sling has two straps with movable Velcro buckle closure for more positive positioning. Swathe can be positioned over sling to anchor arm to body around chest, through sling to anchor arm to body around waist.

■ Indication

For bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



FOAM CLAVICLE SPLINT

515



- This brace reduces stress on the clavicle and limits forward and downward movement of the shoulders. Made of breathable polyurethane foam covered with white cotton stockinet. Two straps with Velcro closure attached through triangular metal piece. Current design features frontal hook and loop closures. 5 sizes

■ Indication

Fractures of the clavicle



RIB BELT

810



- Provides stable support to ribs and torso. Limits the expansion of chest to promote healing. Constructed from Breathable elastic material. Helps to support the ribcage during work or leisure activities.

■ Indication

Post-op, Rib fractures, Sternum Fractures



628-629

BELLY BELT

- Produced from a special rubber to apply pressure on the stomach of newly born babies. It has a soft rounded cushion placed on the belly area. One size. Adjustable. It is a hernia cushion with a soft choke produced from a special rubber for the children. Adjustable. Available in three models as Right-left and a double sided. Only one size.

Indication

- Model 628 belly ruptures
- Model 629 For hernia in the groin area



634-636

INGUINAL RUPTURE BELT

- Cushion that opposes the hernia area is made of a rubber mixture soft material. The part that rounds the waist and passes between the legs is produced of a special elastic rubber band and fits all the adult sizes. Model 634, is produced both as right and as left model for one sided hernias, and the model 636 is produced for double sided hernias in one size only.

Indication

- Waist hernias (Inguinal hernia)



632

INGUINAL HERNIA CORSET

- By means of the two pockets both on the right and left sides it can be used for both left and right hernias and for the double sided hernias as well. Besides the pockets may be filled with require amount of soft material and that allows required pressure on the patients hernia. It is recommended to be worn in laying position.

Indication

- Hernias in the waist area (inguinal hernia)



638

SUSPENSORY

- The part that is opposite the scrotum is made of a special cotton fabric and a soft elastic belt is added for holding the waist.

Indication

- It is used in the non-responsive hernias or injuries in the testicles and in all branches of sports for supplying more freedom of movement and to control the upward pull of the pouch



WOOL KNEE WARMER

420



- Wool provides comfort and warmth with natural fibers and relieves the pain.

Indication

Rheumatism, Pains around knee, when the knee needs stay warm



WOOL ELASTIC BELT 27CM

611



- This Corset is made from Wool. Provides comfort and protects the body temperature of the area which it include, Helps relax the muscles. Tubular and seamless. Long-lasting and anti-rolling. Non-allergic. Comfortable to wear and unnoticeable under clothes. Provides light support and protects lower back, abdominal and hip areas from cold and humidity. Wool provides natural warmth and makes it very comfortable to wear on a daily basis.

Indication

For prevention and treatment of colitis, arthritis, muscle pain, minor injuries. For postoperative and postpartum recuperation. For protection of lower back and abdominal areas during severe weather conditions



WOOL ELASTIC BELT (CURVED) 32CM

612



- This model is five cm longer than the model 611 and makes the same effects. For ease of use its curved on front side.

Indication

For prevention and treatment of colitis, arthritis, muscle pain, minor injuries. For postoperative and postpartum recuperation. For protection of lower back and abdominal areas during severe weather conditions



WOOL ELASTIC BELT (LEG TYPE) 32CM

613



- Same effects with model 611 and 612. It's produced mainly for gynecologic cases also used for sacral, lumbar, abdomen, and perineal.

Indication

For prevention and treatment of colitis, arthritis, muscle pain, minor injuries. For postoperative and postpartum recuperation. For protection of lower back and abdominal areas during severe weather conditions



650-652

ELASTIC ABDOMEN CORSET

- Made from double layer lycra with satin lamination on front side. Wraps the stomach and supports the area.

Indication

For women give birth many times, supporting the waist area, post-op



654-656

ELASTIC ABDOMEN CORSET

- Made from double layer lycra with satin lamination on front side. Wraps the stomach and supports the area.

Indication

For women give birth many times, supporting the waist area



658-660

ELASTIC ABDOMEN CORSET

- Made from double layer lycra with satin lamination on front side. Wraps the stomach and supports the area.

Indication

After pregnancy as a postpartum abdominal support (especially after C-Section); After any other surgery in the abdominal area



644

POST-OP ABDOMEN CORSET

- Supports the area after surgery, Also can be used for women after several births.

Indication

After pregnancy as a postpartum abdominal support (especially after C-Section); After any other surgery in the abdominal area



STOMACH CORSET WITH CUSHION

640



- Decreases pressure and provides excellent support to abdomen and waist, Provides "Body Shaping" effect. This model is for weak bodys.

■ Indication

In any situation where abdominal support is needed



STOMACH CORSET

641-642



- Decreases pressure and provides excellent support to abdomen and waist, Provides "Body Shaping" effect. With 3 flexible metal stays and 5 reinforcement straps. This model is for large bodys.

■ Indication

In any situation where abdominal support is needed



KIDNEY-STOMACH CORSET

661-662



- Decreases pressure and provides excellent support to abdomen, waist and lumbo-sacral areas. Provides "Body Shaping" effect. With cushions on the front side for stomach and sideways for kidneys

■ Indication

In any situation where abdominal support is needed



LUMBAL SUPPORT CORSET (double pull)

826 - 827



- Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back. Two additional side pulls are designed for better fitting and tension adjustment. Cotton inside lining helps minimize sweating and reduces allergic reactions. In model 827 there are for metal stays instead of plastic.

■ Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine



826-1

ELASTIC BACK SUPPORT

- Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back.

Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine. Recommended to those, who lift and move heavy objects




828

LUMBOSTAD CORSET

- This lumbar sacral support helps provide support and compression to the lower lumbar region and abdominal area. This lumbar sacral is ideal for shorter or geriatric patients.

Indication

Recommended for rehabilitation and recovery of serious lower back problems caused by accidents and post-operative conditions, Chronic local or pseudo-radicular lumbalgia accompanying degenerative conditions and muscular and ligament inadequacy



821-824

LUMBA SACRAL CORSET

- High quality five metal closures that adjust for the most comfortable and easy fit. The brace features a 3-pull, lateral lace adjustment, two or four removable, contoured, metal paraspinal stays, anterior hook closure, and mesh gussets at the rib and perineal area.

Indication

Chronic local or pseudo-radicular lumbalgia accompanying degenerative conditions and muscular and ligament inadequacy



831-834

DORSO LUMBAR CORSET

- High quality five metal closures that adjust for the most comfortable and easy fit. It features a 3-pull, side lace adjustment, two removable contoured metal spinal supports, flexible front, side and back stays, and removable shoulder straps.

Indication

Recommended for compression fractures of the thoracic spine, postoperative convalescence, arthritis, osteoporosis, or postural defects such as scoliosis



HYPEREXTENSION BRACE

840



- Designed to act as three point principle by exerting pressure against sternum and pubis on anterior aspect and the Lumbar Region on the Posterior Aspect. This fixes the Lumbar vertebral column in a hyperlordotic position. Easy to apply, lightweight and comfortable.


Indication

Compression Fractures of the lower thoracic and upper lumbar vertebral column. Undisplaced stable fractures Adolescent kyphosis, Any problem of the Thoraco Lumbar Junction which needs moderate immobilisation



TAYLOR BRACE

845



- This brace is designed to maintain hyperextension by limiting anterior and posterior motion as well as rotation. Two paraspinous bars and two lateral uprights attached. Abdominal support secured to lateral uprights with corset lacing. Auxiliary straps with slide buckle closure.

Indication

Compression Fractures of the lower thoracic and upper lumbar vertebral column. Undisplaced stable fractures Adolescent kyphosis, Any problem of the Thoraco Lumbar Junction which needs moderate immobilisation



POSTUREX

805



- Can be worn unnoticed under most clothing. Non-elastic white cotton canvas and three flexible stays enable the Posture Corrector to deliver support to the upper back while promoting proper posture. Shoulder Straps keep the support snug and comfortable while helping ensure that the shoulders are held in a neutral position.

Indication

Helps to improve posture



POSTURE CORRECTOR

806 - 806B



- Can be worn unnoticed under most clothing. Designed to help provide good posture. Pull through strap construction offers more support and allows desired strap adjustment. Constructed of durable and breathable elastic material.



- Model 826 B: with two flexible stays for added support.

Indication

Helps to improve posture



801

TROCHANTER BELT

- The Trochanter belt is used for instability of the sacroiliac joints. It features adjustable side pulls over the sacral area for added tension and stability.

Indication

Sacroiliac joint instabilities



807

THORAX TRACTION BELT

- It's a corset made of a special fabric adjustable with Velcro. It's applied to thorax and used to prevent the patient slipping down from the table when using pelvic traction belt.

Indication

Recommended to those, who needs total traction




808-809

PELVIC TRACTION BELT

- 808 Model is pulled down by one strap. 809 model is pulled down from the sides by double straps. The model 809 is seen in the picture. Suitable for use in clinics surgeries and if needed in homes as well.

Indication

Recommended to those, who needs total traction



812K-812B

SKIN TRACTION KIT

- Foot / Skin traction kit is used as a below knee traction. Used to provide traction to limbs after injury & surgery. The skin traction strips provide maximum aeration to the skin while securely holding its centre extremity to the spreader at the foot end. Elasticated long restraint bandage provides firm fixation to all sides of limb.

812K for child
812B for adult

Indication

When below knee traction is needed



❖ CERVICAL TRACTION KIT

900



- This cervical traction kit allows you to enjoy the benefits of head/neck traction therapy at home, and needs nothing more than a door and a chair or stool. Kit includes a sturdy overdoor steel bracket with pulleys, cord, weight bag, and a machine washable head halter .

■ Indication

When cervical traction is needed



❖ CERVICAL COLLAR (CHILD)

910



- For children ages up to 15, Constructed of single piece heavy plastic, perforated for ventilation. Vinyl padded edges, Rear contact closure.

■ Indication

Cervical syndromes, Cervical Disc Disease



❖ CERVICAL COLLAR

911



- Constructed of double piece heavy plastic, Height easily adjustable , perforated for ventilation. Vinyl padded edges, Rear contact closure.

■ Indication

Cervical syndromes, Cervical Disc Disease, For flexion and hyperextension (whiplash) injuries, torticollis and osteoarthritis



❖ CERVICAL COLLAR WITH CHIN SUPPORT

912



- Constructed of heavy plastic, perforated for ventilation. Vinyl padded edges, Rear contact closure. Height adjustable by simple screw adjustment.

■ Indication

Cervical syndromes, Cervical Disc Disease, For flexion and hyperextension (whiplash) injuries, torticollis and osteoarthritis



914

FOAM CERVICAL COLLAR

- Used for mild cervical lesions, or wherever a reminder against making painful neck movements is required. Made of washable polyurethane foam covered with 100% cotton tubular knit stockinette which obscures soiling and enhances durability for effective longer-term wear. Hook-and-pile fastening.

Indication

Mild Cervical syndromes, Cervical Disc Disease



915

NELSON COLLAR

- Moderate control cervical motion. one piece plastazote foam with Velcro closure for easy application and adjustments.

Indication

Cervical syndromes, Cervical Disc Disease, For flexion and hyperextension (whiplash) injuries, torticollis and osteoarthritis



913

PHILEDELPHIA COLLAR

- Closed cell Plastazote, perforated for ventilation, easy to trim if necessary; rigid anterior and posterior uprights. Bilateral hook-and-pile fastening and adjustment. All components are fully radiolucent. Cradles mandible anteriorly, reaches over occiput posteriorly to enhance rotation control.

Indication

Both for mild to moderate soft tissue damage including flexion and hyperextension (whiplash) injuries, torticollis,, osteoarthritis and relief from muscle tension



916

PHILEDELPHIA COLLAR TRACHEOTOMY

- The unique trachea opening allows healthcare professionals the ability to quickly perform carotid pulse monitoring and emergency tracheotomies while providing cervical immobilization.

Indication

For mild to moderate soft tissue damage including flexion and hyperextension (whiplash) injuries, torticollis,, osteoarthritis and relief from muscle tension. Tracheotomy



❖ CERVICAL BRACE

917



- For milder cervical lesions; excellent choice for hospital or accident site application. Rigid, lightweight frame, incorporating anterior and posterior uprights, mandibular and occipital plates - adjustable in height. Closed cell Plastazote interface between frame and skin, minimizes pressure problems; Cool, lightweight, comfortable to wear.

■ Indication

Cervical syndromes, Cervical Disc Disease, For flexion and hyperextension (whiplash) injuries, torticollis and osteoarthritis



❖ WIRE FRAME CERVICAL COLLAR

919



- A well tolerated, lightweight and comfortable cervical orthosis, especially when control of cervical spine flexion is indicated. Tempered wire frame padded with closed cell Plastazote; both height and circumferential contour are easily adjustable by merely bending by hand. Padded occipital piece, easily positioned and maintained by dual hook-and-pile fastening straps on each side. A remarkable combination of effective motion control and a high degree of wearer acceptance.

■ Indication

For moderate cervical lesions including hyperextension (whiplash), osteoarthritis and cases of mild to moderate soft tissue damage



❖ CERVICAL-THORACIC ORTHOSIS (CTO)

918-1



- Sterno-occipital mandibular immobilizer (SOMI) Provides better control than hard collar, designed to apply forces under the chin and occiput to restrict flexion and extension as well as rotation of the head and cervical spine

■ Indication

Prescribed for cervical conditions requiring a high degree of flexion and rotational control of the cervical spine



■ tracheotomy is possible



❖ CTO COLLAR

918-2




- Same as above, without chest support.

■ Indication

Prescribed for cervical conditions requiring a high degree of flexion and rotational control of the cervical spine



■ x-ray lucent



850

AFO

- AFO is indicated to correct foot drop during ambulation in patients with neurological deficit. Made of durable and resilient polypropylene. Fits into standard shoe. Easy hook and loop closures. Can be easily trimmed or molded with heat.

- **Indication**
Drop Foot, Neurological Deficit



813

FOOT TRACTION KIT

- Foot/Skin Traction Kit is used as a below knee traction. The skin traction strips provide maximum aeration to the skin while securely attaching its centre extremity to the spreader.

How To Use:

Apply foam surface to skin upto the knees. Leave loop space at the foot of the leg. Insert spreader inside this loop space. Tie cord on spreader nose and pass the cord through the pulley and then tie the weight.

- **Indication**
When below knee traction is needed



856

MULTI POSITIONAL FOOT SPLINT

- Nightly antbedsores multipositional brace. Equinus adjustment, Neutral position or Talus. Rotational adjustment.



HEEL ANGLE CONTROL

TALUS : 5° - 11°

NEUTRAL: 0°

EQUINUS: 5° - 11° - 18° - 25°

- **Indication**
Post-surgical of hip, knee, heel and foot. Antialgic. Dorsi-flexor splint at 5°. Antidecubitus. Chronic pain at FASCITIS, cavus



206

ELASTIC BANDAGE

- Porous, elastic, flesh colored tape-allows skin to breathe. Stretches over uneven surfaces.
Available in these widths: 6, 8, 10, 12, 15, 20 cm.

- **Indication**
Prevention or Inhibition of Swelling, Soft Tissue Damages





■ **Model: 465 Hip Abduction Orthosis**

The orthosis is light-weight, sanitary, and moisture-resistant. Designed to maintain the hips in an optimal position for development. Vinyl padded edges.

■ **Indication**

Developmental dysplasia of the hips (DDH) in infants



■ **Model: 470 Hip Abduction Splint**

Made from polyethylene, Vinyl padded edges, fastened by 2 belts coming over the shoulder and with another 2 around the waist. It can easily put on and taken off. It's unaffected by urin, washable.

■ **Indication**

Developmental dysplasia of the hips (DDH) in infants



■ **Model: 475 Pavlik Harness**

Pavlik Harness is one way to treat younger babies (0-8 months) with Developmental Hip Dysplasia (DHD). It is very effective in most cases. The goal of DHD treatment is to hold the head of the femur (top of the thigh bone) within the hip socket to allow the hip to develop normally.

■ **Indication**

Developmental dysplasia of the hips (DDH) in infants



■ **Model: 480 Frejka Pillow**

Therapeutic pillow splint used in treating developmental dysplasia of the hip DDH. The pillow is placed between the thighs of infant to correct congenital hip dislocation by maintaining flexion, abduction and external rotation of the hip.

■ **Indication**

Developmental dysplasia of the hips (DDH) in infants



■ **Model: 485 FlexibleForrester Braun**

Made from rigid metal with leather cover. With assistance of spring equipment, allows the first resistance of patient, then accomplish the necessary abduction. Waist parts and the foot parts are foam padded.

■ **Indication**

Developmental dysplasia of the hips (DDH) in infants



■ **Model: 490 Dennis Brown Splint**

With foot plates. Made of quality aluminum Alloy. Anodized bar. Two bar sizes: No1: Between 31cm - 44cm, No2: Between 45cm - 73cm All plates have protractor marked in 5 degree increments and are clear anodized.

■ **Indication**

developmental dysplasia of the hips (DDH) in infants





■ **Model: K-01**
Sizes between 19-23



■ **Model: K-02** ■
Sizes between 19-28

Hand Made
Genuine Leather
Orthopedic Shoes



■ **Model: K-03**
Sizes between 19-28



■ **Model: K-04** ■
Sizes between 19-23



■ **Model: K-05**
Sizes between 19-23



■ **Model: K-06** ■
Sizes between 24-29



■ **Model: 301-505
Metatarsal Pad**
This product is recommended to help relieve forefoot, and neuroma pain. Erüs molded rubber metatarsal pads are adhesive pads that fit into the shoe.




■ **Model: 330-001
Scaphoid Pad**
Erüs scaphoid pads are perfect for giving flexible arches more support. Simply glue to existing inner sole or new device to reduce shells flexibility.



■ **Model: 330-002
Medial / Lateral Heel Pads**
Correct inversion (rolling out) or eversion (rolling in) of the ankles. 4 sizes, 3 or 5 mm thickness.



■ **Model: 330-006
Torqheel**
Used for making rotation 17 0 inside or outside. For children between 0-7 ages.



■ **Model:310, 311, 312 Aluminium Insole**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made from aluminium and covered with leather




■ **Model:314 Rubber Insole**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made from rubber and covered with velvet



■ **Model:315 Plastic Insole**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made of plastic.




■ **Model:320 Cork Insole**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made of cork.



■ **Model:321 Plastic Insole with Leather covered**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made of plastic with genuine leather covered.



■ **Model:323 Cork Insole with leather covered**

For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made of cork with leather covered.



■ **Model:325 Steel Insole with leather covered**


For lifting the medial longitudinal arch and lateral arch; uplifting of the back part of the root, abduction and adduction of the forefoot. Made of steel with leather covered.



■ **Model:326 Heel Cushions**

This heel cushion from Erüs helps to protect your heel against extra stress. They are suitable for continuous wear and reduce stress to the joints and spine even under constant pressure. By unloading the heel, a shock absorbing effect of 40% is accomplished, not only to joints and the spine but also to muscles, ligaments and the Achilles' tendon.






■ **Model: 327100 Silicone Insoles**

Silicone-Insoles are made of 100% elastic Silicone, anatomically shaped for a perfect fit. Built-in cushions provide a shock absorbing effect to the heel longitudinal arch and the metatarsal area, unloading the entire foot. Silicone-Insoles unload the entire forefoot by 40%. Even distribution of pressure and shock reduces stress to the main pressure areas of the foot, eliminating pain and tenderness.



■ **Model: 327200 Silicone Heelcups with center blue spot**

Silicone-Heelcups are made of 100% elastic Silicone, preformed for maximum comfort. They are suitable for continuous wear and reduce stress to the joints and spine even under constant pressure.



■ **Model: 327210 Silicone Heelcups with medial blue spot**

Silicone-Heelcups are made of 100% elastic Silicone, preformed for maximum comfort. They are suitable for continuous wear and reduce stress to the joints and spine even under constant pressure.



■ **Model: 327220 Silicone Heelcups**

Silicone-Heelcups are made of 100% elastic Silicone, preformed for maximum comfort. They are suitable for continuous wear and reduce stress to the joints and spine even under constant pressure.



■ **Model: 328000 Silicone Insoles 3/4**

High Grade Silicone insoles provide therapy and prevention while offering relief and comfort avoiding pressure and shock waves on the heel and metatarsals. The insole ensures effective cushioning of all ligaments and muscle structures. It supports the arch and is ideal for patients suffering from Achilles Tendonitis, Shin Splits, Heel Pain, Bunions, Flat Feet and Knee Pain.



■ **Model: 327900 Silicone Breast Prosthesis**

Medical grade silicone breast prosthesis. Thought for simple mastectomies. Suitable for both sides.





■ **Model: 327400 Toe Spreaders**

A corrective shield between first and second toes, the extra soft Hydrogel toe spreader effectively prevents deviation of the Hallux into a resulting bunion, while avoiding friction.



■ **Model: 327410 Bunion Shield with ring**

The one step-up in the treatment of symptomatic bunions, this extra soft Hydrogel shield presents with a reinforced loop that guarantees correct fitting. It deflects undue rubbing from the shoe, absorbs lateral shock and prevents friction. A strong elastic fibre is woven within the loop for added strength. The particular design of this orthosis calls for extra care when putting it on and taking it off. Remember, treat gently.



■ **Model: 327420 Bunion guard cushion**

The shell shape HYDROGEL daytime bunion protector is a simple and highly effective design that accommodates the symptomatic bunion as if in a cradle, takes the pressure from the shoe away from the affected spot, and minimizes friction while taking up an insignificant amount of room in the shoe.



■ **Model: 327430 Self Adhesive Bunion Cushions**

The perfect safeguard against troublesome bunions, the extra soft Hydrogel bunion cushion brings about instant relief of the affected area, redistributes lateral shoe pressure, and avoids insidious friction. This hypoallergenic cushion presents with a central cut out for easy accommodation of the subluxated first joint.



■ **Model: 327310 Interdigital wedge metatarsal pad**

An extra soft Hydrogel metatarsal pad with a built-in interdigital wedge to prevent displacement of the orthosis on deambulation. Formulated for the overworked transversal arch of the foot, it also shields the metatarsal heads from encroaching callouses and excessive keratin formation. Strongly recommended for persons who stand on the feet for hours on end.




■ **Model: 327320 Metatarsal Pad with reinforced loop**

This extra soft HYDROGEL metatarsal pad is intended to help eradicate plantar callouses and hard skin, reduces the adverse effects of weight bearing, and alleviates painful rubbing and friction on walking. Specifically designed to treat forefoot soreness, and as a preventive device for those who stand long hours on their feet. An elastic fabric, embedded in the digital loop, reinforces it for added durability, although handling of the orthosis with care is always recommended in order to minimize the possibility of tearing.



■ **Model: 327300 Elastic band with silicone metatarsal insert**

This product follows the design of the more traditional elastic metatarsal band to address met heads problems, both biomechanical, and/or lesions. An extra soft Hydrogel disc has been embedded in the plantar aspect of the elastic band to help the painful metatarsals unload and diffuse pressure. It fits snugly any kind of shoe. Package contains one elastic band.



■ **Model: 327500 Digital Rings**

A one-piece configuration, the extra soft Hydrogel digital ring protects the specific toe(s) afflicted with soft corns, sores, or acts simply as a shield against pressure while taking up a minimum amount of space in the shoe. We reiterate our recommendation to handle this product with care when putting it on or taking it off. Package contains two rings






■ **Model: 327510 Digital Protectors**

This variation on the approach to protecting digits consists of a smooth elastic bandage to which an extra soft Hydrogel disc has been added. It serves the same function as the Hydrogel digital ring. Package contains two rings




■ **Model: 327540 Toe Caps**

The seamless Hydrogel toe cap effectively shelters the distal, dorsal and plantar aspects of the toes, including the nail plate, from all kinds of mild trauma. Its reinforced "hood" provides instant relief for lesions located on the very tip of the toe. Package contains one toe cap.



■ **Model: 327600 Adhesive depression pads**

This extra soft Hydrogel depression pad comes attached to a hypoallergenic surgical tape for easy application anywhere on the foot. Of immense practicality, the goal it achieves is elimination of friction, while it helps prevent the formation of new lesions. The disc itself adheres to the skin so that displacement on walking does not occur. The adhesive strapping can be cut or trimmed depending on where it is to be placed. The maximum thickness of the pad measures 3 mm; its diameter, 2.5 cms.



■ **Model: 327620 Self adhesive corn cushions**

The stick on padding of extra soft Hydrogel has been devised to alleviate the painful callous caused by excessively constricting footwear, and to keep similar conditions in check as it has a prophylactic, preventive effect as well. The central cut out on the oval pad takes direct pressure away from the affected spot and acts as a shield. Totally hypoallergenic. Package contains six pads.



■ **Model: 327640 Self Adhesive Heel Protectors**

The Hydrogel heel protector stands out as an innovative strip adhered to the inner contour of the heel of the shoe to eliminate the irritation and rubbing that occurs not only the "break-in" period of newly acquired footwear, but also with tight-fitting shoes. In addition, it helps avoid the uncomfortable "slip off" effect with shoes one size too large or simply too loose.



■ **Model: 327680 Self adhesive plantar pads**

The Hydrogel plantar pad has fundamentally a preventive application for reducing the painful effects of plantar callouses, while discouraging the formation of new keratotic tissue. Its hypoallergenic backing ensures firm adhesion, and a soothing influence on the plantar aspect of the foot.

Package contains two plantar pads.




■ **Model: 327700 Toe Crest with Loop**

A variation on the SILICONE TOE PROP, the toe crest version presents with a reinforced digital loop, and addresses the same kind of "hammer toes" conditions. The reinforcement of the loop consists of an elastic fabric woven within the silicone for added strength and durability. As with all splints that incorporate a loop, we emphasize the need for handling delicately. Available in right and left. Package contains one toe crest.




■ **Model: 327800 Forefoot cover with silicone pad**

This versatile forefoot cover is furnished with a cushioning extra soft Hydrogel pad that translates into instant comfort for overworked metatarsal heads and collapsed transversal arch. A gentle massage caresses the met-heads and tips of your toes, giving lasting relief to sensitive feet. Made of 100% cotton, it comes in an attractive off-white colour. Hand wash and air-dry, then reuse.



neoprene supports

whenever you need support



WEIGHT SET

AL-1000



- Made of neoprene .Has six pockets each containing lead weights. The weight can be increased and decreased on demand.

AL-1000 = 1,5 KG (6X250gr)
AL-1001 = 2,5 KG (6X420 gr)
AL-1002 = 5,0 KG (6x830gr)

■ Indication

It is used in exercises for strengthening the weak muscles in hands and feet, in sport activities and as a helper in physical therapy



WEIGHT SET

AL-1010



- Made of impertex .Contains steel granules forming the weight. Constant in weight.



■ Indication

It is used in exercises for strengthening the weak muscles in hands and feet, in sport activities and as a helper in physical therapy



ANKLE SUPPORT

AL-2118



- Made of neoprene. It reduces the pain and relieves the muscles and tendons by preserving the body heat. Supports the walking and sports activities without restricting the motion.

■ Indication

Reduces the pain and supplies an effective support in the hurting and injuries around the ankle



ANKLE SUPPORT WITH VELCRO

AL-2128



- Made of neoprene. Easily applied to the patient with the help of Velcro bands.

■ Indication

Reduces the pain and supplies an effective support in the hurting and injuries around the ankle



AL-2135

ANKLE GUARD

- Made of neoprene. It has two anatomically appropriate spiral stays that pass through both the inner and outer parts of the ankle area and has an eight-shaped elastic bandage surrounding the foot to apply the ideal required pressure on the ankle.

■ Indication

Used in the medical treatment of weak and middle intensity twist injuries and as a support for restarting the work and sport activities after the twist injuries. It supplies effective compression in hyperlaxity events



AL-2136

ANKLE SUPPORT WITH STRAP

- Made of neoprene and has an eight-shaped elastic bandage surrounding the foot to apply the ideal required pressure on the ankle.

■ Indication

Used in the medical treatment of weak and middle intensity twist injuries and as a support for restarting the work and sport activities after the twist injuries. It supplies effective compression in hyperlaxity events



AL-2140

ANKLE SUPPORT

■ **MALLEOGUARD**

Malleoguard is very comfortable to wear and streamlined to fit in all kinds of walking shoes and sports shoes. The anatomically contoured Malleoguard lies flat against the side of the joint, thus preventing the ankle from coming out of the malleolar mortise (forward movement of the talus) and protecting the relevant capsular ligaments.

The high degree of functional stabilization of the ankle counteracts ankle twisting without restricting freedom of movement.

■ Indication

Early functional treatment for injuries of the lateral malleolar ligaments (and the bifurcate ligaments) Capsular ligament strains, Chronic ligament insufficiency, Postoperative rehabilitation



AL-2146

ANKLE SUPPORT

■ **ACHILLOGUARD**

The Achilloguard has a viscoelastic heel-wedge to raise the heel for relief of Achilles tendon. A special insert on the back of the Achillotrain leaves tendon pressure-free and provides intermittent compression for the reduction of swelling and edema. Additional heel-wedge provided for unaffected foot thus preventing leg length discrepancies.

■ Indication

Achilles tendonitis, Acute Achilles paratendonitis, Achillobursitis, Post-operative treatment after Achilles tendon ruptures, Posterior Achilles tendon bursitis



ANKLE SUPPORT WITH MALLEOL SHELL

AL-2142



■ **MALLEOSTAB**

Gently supports the ankle by providing controlled compression throughout the joint, while relieving pressure over the malleols. a contoured special insert offers gentle, intermittent compression to the soft tissue for the reduction of swelling and edema promotes proprioception, lightweight, breathable knit provides optimum patient compliance.

■ Indication

Early functional treatment for injuries of the lateral malleolar ligaments (and the bifurcate ligaments) Capsular ligament strains, Chronic ligament insufficiency, Postoperative rehabilitation



ANKLE STABILIZATOR

AL-2152



- Ultra lightweight functional ankle brace. Designed to help support and protect unstable or injured ankles. Allows normal dorsal and plantar flexion while controlling inversion and eversion. Foam ankle padding are combined to provide compression and support.

■ Indication

Early functional treatment for injuries of the lateral malleolar ligaments (and the bifurcate ligaments) Capsular ligament strains, Chronic ligament insufficiency, Postoperative rehabilitation



CALF SUPPORT

AL-2180



- The Calf Skin provides the calf with constant compression for both calf and shin area. Excellent in helping reduce shin splints, Bio Skin's light weight and low profile . Focuses compression on the soft tissue next to the tibia (where you need it) and keeps pressure off of the calf (where you don't need it). It can also be worn over the upper portion of the tibia to alleviate symptoms of Anterior Shin Splints.

■ Indication

Shin Splint, Muscle ruptures, Painful conditions



THIGH SUPPORT

AL-2336



- General purpose support compresses and supports soft tissue while retaining therapeutic heat. Helps reduce muscle strains and pulls. Reversible thermal neoprene support with nylon facing on both sides. The thigh Support gives optimal compression to the quads and hamstrings.

■ Indication

Muscle ruptures, painful conditions, strains of the quadriceps or hamstrings



AL-2235

KNEE GUARD

- Easily slips on to provide support to the knee and therapeutic warmth. Soft nylon lining is comfortable against the skin. Compressive neoprene limits swelling. Closed patella design.

Indication

Swelling and edema of the knee, ligament instability, knee instability, recurrent knee pain, arthritis, degenerative joint disease



AL-2252

POST-OP KNEE GUARD

- Perfect for functional or light activities, this brace allows for sufficient flexibility while still protecting and stabilizing the knee. Medial and lateral composite hinges provide support and help prevent hyperextension. Overall brace weight is extremely lightweight. Easy to apply wrap front thigh and calf closures.

Indication

Ideal for post-surgery or mild ligament injury, knee instability, or degenerative joint disease



AL-2254

POST-OP HINGED KNEE

- Guard Against Lateral Impact. Features a solid, composite lateral hinge. Protects the knee from injury due to lateral side impacts. The wide, composite hinge has no exposed metal and is designed with a stop to prevent hyperextension. Wrap style of the calf and thigh straps allows for easy application and removal. The ideal brace for athletic use in sports such as football, soccer, and skiing.

Indication

Ideal for post-surgery or mild ligament injury, knee instability, or degenerative joint disease



AL-2275

PATELLAR TENDON STRAP

- This Model is easy, simple, comfortable and doesn't hinder activity. It is economical. It applies pressure over the patellar tendon. Strap functions dynamically as knee flexes and extends. Easy to apply and adjust. Strap does not impede circulation or restrict mobility.

Indication

Indicated for chondromalacia



PATELLAR STABILIZER

AL-2285A



- This model has a patellar cutout for correct positioning and movable horseshoe buttress--mounted with Velcro for infinite positioning. Two lateral and two medial spiral stays. Opposing loop-and-lock straps above and below the patellar opening.

■ Indication

For post operative rehabilitation, patellar tendonitis and dislocations, anterior/posterior cruciate instability and lateral/medial collateral ligament instability



HINGED KNEE BRACE

AL-2285B



- This open-patella knee brace is constructed of sports neoprene for stabilizing compression and therapeutic warmth. With steel support hinges to help prevent hyper-extension and provide side-to-side support. Built-in superior horseshoe support controls patellar malalignment. Two adjustable straps keep the brace in place.

■ Indication

Ideal for knee ligament tears, sprains and chronic knee instability problems



PATELLAR KNEE GUARD

AL-2285C



- Easily slips on to provide support and therapeutic warmth to the knee. Soft nylon lining is comfortable against the skin. Compressive neoprene helps limit swelling. Open patella design.

■ Indication

Swelling and edema of the knee, ligament instability, knee instability, recurrent knee pain, arthritis, degenerative joint disease



NEOPRENE KNEE STABILIZER

AL-2290



- Great Fit, Firm Stability. Built-in spring light stays on either side of this support provide side-to-side stabilization of the knee joint. Removable horseshoe insert increases the support around the kneecap and can be placed exactly where needed. Four-way stretch sports neoprene gives proper fit and therapeutic warmth. Two adjustable straps keep the stabilizer in place.

■ Indication

Ideal for knee ligament tears, sprains and chronic knee instability problems



AL-2291

HINGED KNEE WITH DIAGONAL STRAP

- Easy wrap design for easy application and removal. Steel hinges provide medial/lateral support and a removable rubber donut can be cut to a superior or inferior horseshoe to stabilize the kneecap. Diagonal Straps can be easily adjusted for the optimum pressure on the ligaments. Constructed of sports neoprene.

Indication

Ideal for knee ligament tears, sprains and chronic knee instability problems



AL-2295

R.O.M KNEE BRACE

- Rapid application - universal size fits all. Toolless joint for quick and accurate R-O-M control in 10° flexion/extension increments. Simple length adjustment via slotted hinge arms. Aluminum hinge arms may be contoured at condyles for optimum fit.

Indication

For post operative rehabilitation, patellar tendonitis and dislocations, anterior/posterior cruciate instability and lateral/medial collateral ligament instability



AL-2296

LIGHTWEIGHT POLYCENTRIC HINGED KNEE BRACE

- This hinged knee sleeve features lightweight but durable dual-axis hinges that allow for a full range of motion. Built-in stops in the hinges prevent hyperextension, and contoured shells provide additional protection against lateral impacts during activity. Open patella design with four adjustable support straps. The low profile, slip-on style is easy to apply and is comfortable to wear.

Indication

For post operative rehabilitation, patellar tendonitis and dislocations, anterior/posterior cruciate instability and lateral/medial collateral ligament instability



AL-2299

KNEE IMMOBILIZER

- Provides immobilization of the knee after surgery or injury. It has one medial and one lateral stay, plus three posterior contoured stays. Medial and lateral stays are movable to accommodate different knee circumferences. The immobilizer straps onto the knee with loop lock closures.

Indication

When immobilization needed after surgery or injury



RIB BELT

AL-2454



■ **COSTAFIX**

Provides stable support to ribs and torso. Limits the expansion of chest to promote healing. Constructed from Breathable elastic material. Helps to support the ribcage during work or leisure activities.

■ Indication

Post-op, Rib fractures, Sternum Fractures



ABDOMEN CORSET

AL-2484



■ **ABDOMITEX**

Provides stable support to ribs and torso and abdominal area. Limits the expansion of chest to promote healing. Constructed from Breathable elastic material. Helps to support the ribcage during work or leisure activities.

■ Indication

Post-op, Rib fractures, Sternum Fractures, when abdominal support is needed



NEOPRENE BACK SUPPORT

AL-2480



■ Made of neoprene. Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back.

■ Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine. Recommended to those, who lift and move heavy objects



ABDOMINAL BACK SUPPORT

AL-2485



■ Constructed from Breathable elastic material. Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back.

■ Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine. Recommended to those, who lift and move heavy objects



AL-2486

LUMBO SACRAL CORSET

LUMBOLIGHT

Constructed from Breathable elastic material. Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back. Back height is 26 cm, covers the area from sacrum to D12

Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine. Recommended to those, who lift and move heavy objects



AL-2490

LUMBA SACRAL CORSET (tall level)

LUMBOTEX

Constructed from Breathable elastic material. Provides excellent support, warmth and comfort to the abdominal and lumbo-sacral areas. Keeps your lower back warm, protects from injuries and helps reduce lower back pain. Four flexible plastic stays provide increased stability and support to the lower back. Back height is 32 cm, covers the area from sacrum to D9.

Indication

Recommended to those, who suffer from osteochondrosis, lumbago or have ever experienced lower back pain, who have undergone surgery on the lower section of the spine. Recommended to those, who lift and move heavy objects



AL-2495

DORSO LUMBAR CORSET

DORSOTEX

Constructed from Breathable elastic material. Provides excellent support, warmth and comfort to the abdominal and dorso-lumbar areas. There are four flexible stays, two of them from sacrum to D9, and two of them from sacrum to D4, and there are two rigid stays from sacrum to D4. Elastic shoulder straps provides extension support and postural correction.

Indication

Recommended for compression fractures of the thoracic spine, postoperative convalescence, arthritis, osteoporosis, or postural defects such as scoliosis



AL-2498

POSTURAL EXTENSION CORSET

POSTUREFIX

Constructed from Breathable elastic material. Can be worn unnoticed under most clothing. Two flexible stays enable the Posturefix to deliver support to the upper back while promoting proper posture. Shoulder straps keep the support snug and comfortable while helping ensure that the shoulders are held in a neutral position.

Indication

Helps to improve posture



❖ SACROBELT

AL-2460



■ **SACROBELT**

The Sacrobelt is designed to compress and support the sacroiliac joints, thereby relieving stress and instability at these weight bearing structures. Just as importantly, it is not so tight that it is restrictive. Excess or insufficient motion can adversely affect the entire musculoskeletal system. The Sacrobelt, by providing the correct balance of resistance and resilience, re-establishes the joint's normal motion.


■ Indication

Sacroiliac joint instabilities



❖ SACROBELT THERMO

AL-2465



■ **SACROBELT thermo**

This corset is designed to compress and support the sacroiliac joints, thereby relieving stress and instability at these weight bearing structures.

■ Indication

Sacroiliac joint instabilities



❖ MATERNITY SUPPORT

AL -2475



■ This corset supports in a natural way. It lifts and supports the abdomen and transfers the weight evenly and comfortably to spinal system. It fastens quickly with simple Velcro closures and is designed to grow along with you so you can wear it right through to delivery.

■ Indication

Pregnancy



❖ THERMAL WEAR

AL-2442



■ These neoprene shorts provide warmth and support to the back, stomach, hip flexer and thighs and are ideal for use with hydrotherapy or exercise programs.

■ Indication

Helps to loose weight, helps physiotherapy exercises



AL-2575

WRIST SUPPORT

- Made of Neoprene. Help maintain normal hand activity after sprains and with tendinitis. Used to help prevent strain and injury during work or athletics. Easy off - Easy on velcro closure.

- Indication**
Sprains, Tendinitis



AL-2584

THUMBOSTAB

- THUMBOSTAB thermo**
Helps support and limit the motion of the MP joint of the thumb. Ideal for use on soft tissue injuries, ligament strains, Gamekeeper's thumb, osteoarthritis, and degenerative joint disease. Allows for full finger dexterity. Thermal lining provides warmth to the thumb area and joint..

- Indication**
Gamekeeper's Thumb, CMC and MC injuries



AL-2585

WRIST SPLINT

- Made of Neoprene. It's used when the wrist is strained and twisted or when the wrist needs to have a rest. Designed to stabilize the wrist in the neutral position. This splint stabilizes the wrist while still allowing full movement of the fingers. Features a removable palmar stay for extra support and immobilization. It can be used on volar or dorsal sides, this brace have two stay bed, it can be used for left or right.

- Indication**
Rheumatic and traumatic arthritis, tendinitis and carpal tunnel, post op



AL-2586

MANUTHUMBSTAB

- MANUTHUMBSTAB thermo**
Made of Neoprene. Recommended when immobilization of the wrist and/or thumb is needed. Provides intimate support, immobilization and thumb abduction.

- Indication**
Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, Gamekeeper's Thumb, CMC and MC injuries



THUMBOSTAB

AL-2594



■ **THUMBOSTAB** airtex

Helps support and limit the motion of the MP joint of the thumb. Ideal for use on soft tissue injuries, ligament strains, Gamekeeper's thumb, osteoarthritis, and degenerative joint disease. Allows for full finger dexterity. Thermal lining provides warmth to the thumb area and joint..

■ Indication

Gamekeeper's Thumb, CMC and MC injuries



MANUSTAB

AL-2590



■ **MANUSTAB** airtex

Constructed from Breathable elastic Airtex Neoprene. It's used when the wrist is strained and twisted or when the wrist needs to have a rest. Designed to stabilize the wrist in the neutral position. This splint stabilizes the wrist while still allowing full movement of the fingers. Features a removable palmar stay for extra support and immobilization.

■ Indication

Rheumatic and traumatic arthritis, tendinitis and carpal tunnel, post op



MANUTHUMBSTAB

AL-2596



■ **MANUTHUMBSTAB** airtex

Constructed from Breathable elastic Airtex Neoprene, recommended when immobilization of the wrist and/or thumb is needed. Provides intimate support, immobilization and thumb abduction. Designed to stabilize the wrist in the neutral position. Features a removable palmar stay for extra support and immobilization.

■ Indication

Rheumatic and traumatic arthritis, tendonitis and carpal tunnel, post op, Gamekeeper's Thumb, CMC and MC injuries



ELBOW SUPPORT

2552A



■ This Elbow Support is designed to provide superior compression with a neoprene construction. This comfortable brace aids in the prevention and maintenance of elbow injuries.

■ Indication

Elbow injuries



AL-2552B

ELBOW SUPPORT REINFORCED

- This Elbow Support is designed to provide superior compression with a neoprene construction. It has a 5 cm wide elastic strap on the lateral epicondylitis area. This comfortable brace aids in the prevention and maintenance of elbow injuries.

■ Indication

Elbow injuries, Epicondylitis



AL-2554

EPICONDYLITIS BANDAGE

- This armband applies focused compression to the tendon without restricting circulation. Eliminates stress and helps reduce inflammation and pain. Ideal for wear during athletic activity, work, or at rest.

■ Indication

Epicondylitis, tennis elbow and muscle strains



AL-2556

EPISPORT

■ **EPISPORT**

Helps relieve pain from muscular overuse caused by work, sports or injury. Provides targeted compression to the long tendons of the forearm allowing them to remain parallel and pain free.

■ Indication

Epicondylitis, tennis elbow and muscle strains



AL-2626

SHOULDER SUPPORT

- Made of neoprene. Provides warmth, compression and support to shoulder and aids healing of the shoulder. Velcro closure for pulling and off easily.

■ Indication

Sprains, strains, bursitis and tendonitis



FULL SHOULDER SUPPORT

AL-2627



- Made of neoprene. Provides warmth, compression and support to shoulder and aids healing of the shoulder. Velco closure for pulling and off easily.

■ Indication

Sprains, strains, bursitis and tendonitis



ELASTIC SHOULDER IMMOBILIZER

AL-2632



- Immobilizes shoulder comfortably after separation, dislocation, strains and other shoulder injuries. Forearm cuff, humeral cuff, and shoulder strap may be placed in desired position. Shoulder strap with velcro closure.

■ Indication

Provides immobilization for bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



ACROMIOCLAVICULAR SPLINT

AL-2642



- Padded shoulder and chest straps ,Foam padding in pouch for greater patient comfort, Hook and loop closure helps secure arm position ,Buckle adjustments, Easy application ,Fits left or right .

■ Indication

Provides immobilization for bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



CLAVICLE SUPPORT

AL-2644



- This brace reduces stress on the clavicle and limits forward and downward movement of the shoulders. High stability provided by non-slip strap material. Quick and easy to fit. Encourages the correct alignment of the fractured bone by pulling the shoulders backwards

■ Indication

collarbone dislocations, clavicle fractures



AL-2682

VELPOTEX

■ VELPOTEX

Immobilizes shoulder comfortably after separation, dislocation, strains and other shoulder injuries. Hook and loop closure helps secure arm position.

■ Indication

Provides immobilization for bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



AL-2690

FOAM ARM SLING

- Superior arm sling that evenly distributes weight, while holding arm securely in place. Made of durable foam padded cotton with hook and loop closure for easy adjustment and maximum comfort.

■ Indication

Sprained, broken or surgically operated arm while it recuperates



AL-2692

ArmSling 30° Abduction

- Can be used after rotator cuff repair .Helps keep the shoulder in slight 30° abduction.Used as a transition between immobilization with a shoulder abduction splint and full function without a splint .

■ Indication

For bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



AL-2694

Arm Sling 45° - 60° Abduction

- Can be used after rotator cuff repair .Helps keep the shoulder in slight 45° or 60° abduction. Provides post-injury or post-surgical support of arm and shoulder.

■ Indication

For bankhart repairs, rotator cuff repairs, capsular shifts and other shoulder instabilities



Hot - Cold Compress

THERMO JEL

AL-1100




- It preserves the heat along 30 minutes. Should be stored in freezer 1 hour for cold application or ten minutes in hot water for heat therapy. It can be easily adjusted with elastic velcro straps, and can be fixed in the requisite area of the body.

Indication

Cold Therapy: Various disorders

Heat Therapy: Rheumatism and muscle pains



THERMO GEL

AL-1120



- It preserves the heat along 30 minutes. Should be stored in freezer 1 hour for cold application or ten minutes in hot water for heat therapy.

Indication

Cold Therapy: Various disorders

Heat Therapy: Rheumatism and muscle pains



First Aid Splints

AIR SPLINTS

AL-1200

AL-1200



Full Set (8 Unit)

- Our inflatable air splint has a long list of strengths. It is easy to apply and remove. It has a push-pull valve for hassle-free inflation. It's x-ray transparent. And it comes with an available zipper closure. This air splint provides strong and reliable support for fractured or broken bones.



AL-1210	Foot & ankle	25,4cm
AL-1220	Foot & ankle	38,1cm
AL-1230	Half Leg	63,5cm
AL-1240	Full Leg	81,2cm



AL-1250	Hand or Wrist	25,4cm
AL-1260	Hand or Wrist	38,1cm
AL-1270	Half Arm	63,5cm
AL-1280	Full Arm	81,2cm

ALTEX[®]
Medical Stockings



ANKLE SUPPORT

1050



■ This Ankle Support is constructed of heavy-duty elastic. It's ideal for warm weather climates or patients allergic to neoprene, for mild soft tissue injuries and arthritis. Provides longlasting and effective compression.

■ Indication

Tendon injuries, Ankle sprains



KNEE SUPPORT

1055



■ This Knee Support is constructed of heavy-duty elastic. It's ideal for warm weather climates or patients allergic to neoprene. Provides longlasting and effective compression.

■ Indication

For mild soft tissue injuries and arthritis



BELOW KNEE

1010



■ Firm support from ankle to calf, heeled, closed toed. Flesh coloured. Durable and Effective. Provides longlasting and effective compression.

■ Indication

Varicosis marked by a tendency to oedema. Post-traumatic swelling. Following superficial thrombophlebitis. Following sclerosis and surgical treatment of varicose veins. Very severe pre-natal varicosis. Following treatment for light ulcers



1020

THIGH HIGH

- Thigh-length stocking, fine weave. Firm support from ankle to thigh, heeled, closed toed. Flesh coloured. Durable and Effective. Provides longlasting compression.

Indication

Varicosis marked by a tendency to oedema. Post-traumatic swelling. Following superficial thrombophlebitis. Following sclerosis and surgical treatment of varicose veins. Very severe pre-natal varicosis. Following treatment for light ulcers




1030

PANTY


- Designed to prevent venous disorders of the leg such as varicose veins, spider veins and ankle/foot swelling. Provides moderate graduated compression to increase circulation in the legs.

Indication

Varicosis marked by a tendency to oedema. Post-traumatic swelling. Following superficial thrombophlebitis. Following sclerosis and surgical treatment of varicose veins. Very severe pre-natal varicosis. Following treatment for light ulcers



izidor[®]
connection to life



Model: 710 Wooden Crutches

Manufactured from wooden laminate of the highest quality with a foam comfy under arm pad, a comfy hand grip and a grey non-slip rubber tip.



Model 712 Aluminium Crutch

Includes tips, grips, underarm with soft foam rubber. Adjustable height.



Model: 718 Cane

Hand walking stick made of aluminium. Adjustable height. Includes tips, spring, button adjustment.



Model: 714 Forearm Crutch

Produced of aluminium. Adjustable height. Includes tips, spring button adjustment.



Model: 722 Fixed Walker

Manufactured from chrome plated steel tubing. fitted with a moulded plastic handgrip and grey non slip rubber tips. fixed height.



Model: 722U With adjustable height

Manufactured from chrome plated steel tubing. Fitted with a moulded plastic handgrip and grey non slip rubber tips. Adjustable height.



Model: 720 Triple Cane Model: 721 Quad Cane

Manufactured from chrome plated Steel tubing. They are all fitted with grey comfy moulded handles and grey non slip rubber tips. All models are suitable for left and right hand use.



Model: 722KU Walker

Adjustable aluminium walking frame which folds easily for storage or transportation.



Model: 723 Walker

Aluminium adjustable reciprocal walking frame which easily folds flat for storage or transportation.





■ **Model: 728 Tri-Wheel Walker**

Chromium-plated steel frame with solid wheels and 'loop' type locking hand brakes. Folds flat for easy storage and transportation.



■ **Model: 730 Rollator**

A rigid frame walker fitted with solid front wheels and two rear rubber tips. Height adjustable.



■ **Model: 930-1 Commode Chair**

Steel tubing with plastic handgrips, Detachable pan with lid, Adjustable height between 43-53 cm, Hygienic in use and easy to clean since it's made of aluminum and plastic



■ **Model: 335 Peroneus Brace**

This brace is used as a support for the feet which is paralysed. In paralytic cases not using peroneus. Is used as a support for paralysed foot. It's adaptable to any kind of shoes and as it is seen in the picture it can be removed and Re-Fitted any time and also can be used with various shoes by having spare connections. 5 sizes available.

Indication:
Drop Foot



■ **Model: 725 Walking Heel Single**

■ **Model: 726 Walking Heel Double**

Made of rubber, highly durable and resilient. Wide platform design for additional walking surface and increased stability and also allows rocking action. Medial longitudinal arch support for better stability. Used for both right and left.



■ **Model: 727 Cast Shoe**

Upper side made of strong canvas duck. Straps of canvas duck have hook and loop closures for a more secure fit. The rocker sole aids in the mobility of the patient and the open toe and heel allows for more versatility in sizing.



■ **Model: 724 Crutch Tips**

Extra wide base for providing maximum stability, fitted with galvanized steel plates to evenly spread pressure across the ferrule prolonging life.

WHEELCHAIR











920K

Color Options



- This model combines the simplicity and style. Offset permanent arms offer full seat width with a permanent arm. Pneumatic Tires. Fixed footrest. Folds flat for easy storage and transportation.



									
65 cm	42 cm	33 cm	42 cm	106 cm	52 cm	96 cm	62 cm	18 cm	18,5


WHEELCHAIR











920KH

Color Options



- If it's quality and durability you're looking for and a moderate price then this low maintenance wheelchair also has a full complement of options and accessories. Its contemporary design and high-performance double cross frame make the 920 KH a definitively superior value. Detachable armrest, detachable footrest. Pneumatic tires.



									
65 cm	42 cm	31 cm	42 cm	106 cm	51 cm	92 cm	62 cm	18 cm, 20 cm	21,5

920KHA

WHEELCHAIR

- This model has all the features of 920 KH plus swinging footrest.




Color Options

65 cm	42 cm	31 cm	42 cm	106 cm	51 cm	92 cm	62 cm	18 cm, 20 cm	21,5

920KHJ

WHEELCHAIR

- This model is for children. It has all the features of 920 KH.



Color Options

62 cm	38 cm	28,5 cm	38 cm	102 cm	51 cm	92 cm	62 cm	18 cm	20,5

WHEELCHAIR RECLINING BACK

920KY

Color Options



- The adjustable-angle back easily changes from 0 degrees to 80 degree. There are three different angle position. You may fix the first position at 12°, second position at 45°, third position at 80° Double Cross, Detachable armrest, detachable footrest, Folds flat for easy storage and transportation.



66 cm	42 cm	34 cm	42 cm	120 cm	52 cm	103 cm	62 cm	18 cm	24,5

WHEELCHAIR

920KO

Color Options



- Its ideal for indoor use, Standard WC, brakes, folds flat for easy storage and transportation.

47 cm	43 cm	15 cm	40 cm	80 cm	50 cm	90 cm	10 cm	13,5

920CC

WHEELCHAIR

- For patients with weak muscles or difficulty with one arm, it can be moved with one arm and directed with other arm. Folds flat for easy storage and transportation.



Color Options

74 cm	42 cm	42 cm	42 cm	106 cm	51 cm	96 cm	62 cm	18 cm	22,5

930

COMMODE CHAIR WITH WHEEL

- Constructed of steel chrome plated tubing. The plastic commode seat is removable and cleanable and the grey plastic pail slides out from underneath the seat from behind. The vinyl padded seat cover goes on top of the commode and converts the commode chair into a transfer chair, four swivel casters with rear caster brakes.



Color Options

54 cm	42 cm	38 cm	51 cm	49 cm	98 cm	10 cm	8,5



Model	Small	Medium	Large	XLarge	XXLarge
301-001	18-19	20-21	22-23	24-25	26
304-021	18-20	21-23	22-23	24-25	26



Model	Small	Medium	Large	XLarge	XXLarge
301-007	18-20	21-23	23,5-25	26-27,5	28,5-30
301-077	27-29	30-34	35-38	-	-



Model	XXsmall	Xsmall	Small	Medium	Large	XLarge	XXLarge
301-006	40-43	43-46	46-49	49-52	52-55	55-58	58-61




Model	Small	Medium	Large	XLarge	XXLarge
508	<3	3-4	4-6	6-7	7-8



Model	Small	Medium	Large	XLarge	XXLarge
504	-	12-16	17-21	-	-
504-001	14-16	16-18	18-20	20-22	-
505	13-14	14-16	16-18	18-20	20-22
505-001	13-14	14-16	16-18	18-20	20-22
506	13-14	14-16	16-18	18-20	20-22
520	13-14	14-16	16-18	18-20	20-22
521	13-16	16-19	19-22	-	-
522	13-16	16-19	19-22	-	-
523	13-16	16-19	19-22	-	-
524	13-16	16-19	19-22	-	-



Model	Small	Medium	Large	XLarge	XXLarge
510	30-36	36-42	42-48	-	-
512	30-36	36-42	42-48	-	-
513	30-36	36-42	42-48	-	-
514	30-36	36-42	42-48	-	-
518	30-36	36-42	42-48	-	-



Model	Small	Medium	Large	XLarge	XXLarge
515	60-75	75-90	90-105	105-115	115-125
807	60-80	80-100	100-120	-	-
918-1	50-70	70-90	90-110	-	-



Model	Small	Medium	Large	XLarge	XXLarge
420	30-33	33-36	36-39	39-42	42-45






Model	1	2	3	4	5	6	7	8
465	14-18	18-22	22-26	26-30	-	-	-	-
470	14-15	16-17	18-19	20-21	22-23	24-25	26-27	28-29
480	14-15	16-17	18-19	20-21	22-23	24-25	26-27	28-29



Model	Small	Medium	Large	XLarge	XXLarge
632	50-70	70-90	90-110	-	-
638	50-60	60-70	70-80	80-90	90-100
801	32-35	35-38	38-42	42-46	46-50
808	28-34	34-40	40-45	-	-
809	28-34	34-40	40-45	-	-



Model	Small-Medium	Large-XLarge
840	Sternum-Pubis → 40-47	48-60
	Gobek Çevresi → 75-90	80-110



Model	Small	Medium	Large	XLarge	XXLarge
611	50-70	70-100	90-110	-	-
612	50-70	70-100	90-110	-	-
613	50-70	70-100	90-110	-	-
640	55-60	60-65	65-70	70-75	75-80
641	60-66	66-72	72-78	78-84	84-90
642	62-69	69-76	76-83	83-90	90-100
644	60-66	66-72	72-78	78-84	84-90
661	60-67	67-74	74-81	81-88	88-95
662	62-69	69-76	76-83	83-90	90-100
805	56-64	64-72	72-80	80-90	90-100
806	56-64	64-72	72-80	80-90	90-100
806B	56-64	64-72	72-80	80-90	90-100
810	75-85	85-95	95-105	105-115	115-125
826	65-75	75-85	85-95	95-105	105-115
827	65-75	75-85	85-95	95-105	105-115
826-1	65-75	75-85	85-95	95-105	105-115
828	50-70	60-80	70-90	80-100	90-110



Model	Small	Medium	Large	XLarge	XXLarge
821	Waist cir. 50-70	60-80	70-90	80-100	90-110
831	Hip cir. 70-90	80-100	90-110	100-120	110-130
	Front height 25	26	27	28	29
	Back height 39 / 55	39 / 55	39 / 55	39 / 55	39 / 55
822	Waist cir. 60-80	70-90	80-100	90-110	100-120
832	Hip cir. 70-90	80-100	90-110	100-120	110-130
	Front height 25	26	27	28	29
	Back height 39 / 55	39 / 55	39 / 55	39 / 55	39 / 55
823	Waist cir. 50-70	60-80	70-90	80-100	90-110
833	Hip cir. 70-90	80-100	90-110	100-120	110-130
	Front height 22	23	24	25	26
	Back height 34 / 45	34 / 45	34 / 45	34 / 45	34 / 45
824	Waist cir. 60-80	70-90	80-100	90-110	100-120
834	Hip cir. 70-90	80-110	90-110	100-120	110-130
	Front height 22	23	24	25	26
	Back height 34 / 45	34 / 45	34 / 45	34 / 45	34 / 45

821 / 831 : Tall Woman
822 / 832 : Tall Man
823 / 833 : Short Woman
824 / 834 : Short Man

■ MEASURING METHOD
Pelvic line should be proper and tight measured. Pay attention to patients sex, man or woman, short or tall.



Small

Normal Woman Neck




Medium

Thin-Tall
Man-Woman Neck



Large

Normal Man Neck



XLarge

Weight-Short
Man Neck




XXLarge

Weight-Tall
Man Neck

Cervical Collar
measurements
should be fit
patients neck



Model	Small	Medium	Large	XLarge	XXLarge
AL-2118	19-22	22-25	25-28	28-31	31-34
AL-2128	19-22	22-25	25-28	28-31	31-34
AL-2135	19-22	22-25	25-28	28-31	31-34
AL-2136	19-22	22-25	25-28	28-31	31-34
AL-2140	19-22	22-25	25-28	28-31	31-34
AL-2142	19-22	22-25	25-28	28-31	31-34
AL-2146	19-22	22-25	25-28	28-31	31-34
AL-2152	19-22	22-25	25-28	28-31	31-34



Model	Small	Medium	Large	XLarge	XXLarge
AL-2180	33-35	35-37	37-39	39-41	41-43




Model	Small	Medium	Large	XLarge
AL-2296	39-47	47-53	53-59	59+
AL-2336	49-51	51-54	54-57	57-60



Model	Small	Medium	Large	XLarge	XXLarge
AL-2235	30-33	33-36	36-39	39-42	42-45
AL-2252	30-33	33-36	36-39	39-42	42-45
AL-2254	30-33	33-36	36-39	39-42	42-45
AL-2275	-	30-36	36-42	-	-
AL-2285A	30-33	33-36	36-39	39-42	42-45
AL-2285B	30-33	33-36	36-39	39-42	42-45
AL-2285C	30-33	33-36	36-39	39-42	42-45
AL-2290	30-33	33-36	36-39	39-42	42-45
AL-2291	30-33	33-36	36-39	39-42	42-45
AL-2295	30-33	33-36	36-39	39-42	42-45
AL-2299	30-33	33-36	36-39	39-42	42-45



Model	Small	Medium	Large	XLarge	XXLarge
AL-2454	75-85	85-95	95-105	105-115	115-125
AL-2484	75-85	85-95	95-105	105-115	115-125
AL-2480	75-85	85-95	95-110	110-125	-
AL-2485	64-72	72-80	80-90	90-100	100-110
AL-2486	64-72	72-80	80-90	90-100	100-110
AL-2490	64-72	72-80	80-90	90-100	100-110
AL-2495	64-72	72-80	80-90	90-100	100-110
AL-2498	64-72	72-80	80-90	90-100	100-110
AL-2475	75-95	95-115	115-130	-	-
AL-2442	65-75	75-85	85-95	95-105	105-120
AL-2682	65-85	85-105	105-125	-	-



Model	Small	Medium	Large	XLarge	XXLarge
AL-2460	86-93	93-100	100-107	107-112	112-119
AL-2475	86-93	93-100	100-107	107-112	112-119






Model	Small	Medium	Large	XLarge	XXLarge
AL-2575	-	14-19	19-24	-	-
AL-2584	13-14	14-16	16-18	18-20	20-22
AL-2585	13-16	16-19	19-22	22-25	-
AL-2586	13-16	16-19	19-22	22-25	-
AL-2590	13-16	16-19	19-22	22-25	-
AL-2594	13-16	16-19	19-22	22-25	-
AL-2596	13-16	16-19	19-22	22-25	-



Model	Small	Medium	Large	XLarge	XXLarge
AL-2552A	23-25	25-28	28-31	31-33	33-35
AL-2552B	23-25	25-28	28-31	31-33	33-35
AL-2554	23-25	25-28	28-31	31-33	-
AL-2556	24-25	26-27	28-30	31-33	-



Model	Small	Medium	Large	XLarge
AL-2626	24-28	20-32	32-35	35-38



Model	Small	Medium	Large	XLarge	XXLarge
AL-2454	75-85	85-95	95-105	105-115	115-125
AL-2627	75-87	87-99	99-111	111-123	-
AL-2632	-	80-100	100-120	-	-
AL-2644	60-75	75-90	90-105	105-115	115-125



Model	Small	Medium	Large
AL-2642	30-36	36-42	42-48
AL-2690	30-36	36-42	42-48
AL-2692	30-36	36-42	42-48
AL-2694	30-36	36-42	42-48



[illegible]

■ Headquarters

Tel: 0.232.264 24 24
Fax: 0.232.254 42 84
E-mail: info@izor.com.tr
Address: 5822 SokakNo: 83 35380 İzmir/Türkiye

■ İstanbul Office

Tel: 0.212.621 49 81-82
Fax: 0.212.621 49 83
E-mail: istanbul@izor.com.tr
Address: Ođuzhan Cad. No: 10 Fındıklı/İSTANBUL

■ Ankara Office

Tel: 0.312.435 12 99
Fax: 0.312.435 29 81
E-mail: ankara@izor.com.tr
Address: Ataç 2. Sokak No: 64/B Kızılay/ANKARA

■ Authorized Dealer